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Laughlin, J. - Microbiology collection device and ISO15189: an appropriate selection - Bacteriology - 2017

ZIEL

Assessment of the performance of a range of swab collection devices with regards to the loss of bacterial viability during significant holding period following sampling. Is subsequent recovery by culture compromised? Transport swab systems investigated:

- MWE Transwab® Amies Gel w/ charcoal (M40)
- Deltalab Amies Gel Swab w/ charcoal
- Copan M40 transport swab Amies Gel w/ charcoal
- Copan non-M40 transport swab Amies Gel w/ charcoal
- TSC Proback Amies Gel w/ charcoal
- Sarstedt Amies Gel w/ charcoal

ERGEBNISSE

The data presented here clearly demonstrate the superior performance of two brands of collection devices, Copan M40-compliant and MWE M40-compliant, in recovery of fastidious organisms such as *S. pneumoniae*, *N. gonorrhoeae* and *H. influenzae*. Overall, all collection devices evaluated appeared to have similar bacterial release at time point zero. Additional concerns have arisen concerning non-viable Gram-negative organisms found in the transport media of Sarstedt and Deltalab collection devices.

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This study illustrates the loss of viability during any significant holding period that compromises recovery. This loss of viability certainly varies significantly between different collection devices. This has significant implications in public health as fastidious organisms or less-robust strains may be underrepresented in the sampled population. In such cases, laboratories may need to assess the suitability of molecular methods for detection of these organisms.

MWE Transwab®:
Neisseria gonorrhoeae held for 24 h and 48 h at 4°C.

Sarstedt device:
Neisseria gonorrhoeae held for 24 h and 48 h at 4°C.

Copan M-40:
Neisseria gonorrhoeae held for 24 h and 48 h at 4°C.

Copan non M-40:
Neisseria gonorrhoeae held for 24 h and 48 h at 4°C.

Laughlin, J. - Efficacy of two commercially available liquid transport collection devices - Poster at the Biomedical Science Congress 2019

ZIEL

Evaluation of the performance of Copan 480C eSwab™ and Medical Wire MW1765 Σ-Transwab® as widely used liquid transport collection devices with regards to their ability to maintain the viability of fastidious strains of clinically significant bacteria and to look at potential overgrowth with *Pseudomonas*.


ERGEBNISSE

The two systems differed in the recovery rates for *Streptococcus pneumoniae* and *Neisseria gonorrhoeae*. Σ-Transwab® demonstrated enhanced recovery rates at the various parameters: *Streptococcus pyogenes*, *Pseudomonas aeruginosa*, *Streptococcus pneumoniae*, *Neisseria gonorrhoeae*, and *Staphylococcus aureus*.

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Copan eSwab™ suffered from significant overgrowth of *Pseudomonas aeruginosa* as a factor of time. There was complete overgrowth in a mixed culture with *Streptococcus pneumoniae* such that colonies could not be detected due to the overgrowth by *Pseudomonas*.

Bacteria and swab type	Temperature	Bacterial recovery (CFU) average of 3			Compliant
		0 hr.	24hr	48hr	
<i>Pseudomonas aeruginosa</i> ATCC 27853					
Sigma Transwab® MW1765	RT	>250	>230	>250	Yes
	4°C	>250	>210	>180	
Copan Diagnostics 480C eSwab®	RT	>250	>250	>250	Overgrowth?
	4°C	>250	>250	>250	

Pseudomonas aeruginosa RT for 24hrs

Laughlin, J. - Molecular efficacy of two commercially available liquid transport collection devices for downstream testing on the GeneXpert® Platform and the Roche FLOW® system - Poster at the European Meeting on Molecular Diagnostics 11th Edition 2019

ZIEL

Evaluation of the performance of Copan eSwab™, Medical Wire Σ-Transwab® Foam tip and Σ-Transwab® PurFlock® tip for the transport and maintenance of Influenza A H3N2 Virus, Influenza B Virus, Respiratory Syncytial Virus (RSV-A) Type A, *Staphylococcus aureus* ATCC 25923, *Streptococcus pyogenes* ATCC 12344 and *Neisseria gonorrhoeae* ATCC 43069 for molecular testing, and to compare its performance with the gold standard transport medium for viruses: Medical Wire Σ-Virocult®.

ERGEBNISSE

There is a significant difference in molecular performance between eSwab™ and MWE's Σ-Transwab®, Σ-Transwab® PF and Σ-Virocult®. All MWE devices are able to detect all targets with significantly lower CT values, meaning MWE devices have greater sensitivity and will give lower limits of detection on molecular platforms.

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The Σ-Transwab® Foam tip and Σ-Transwab® PurFlock® tip are a truly open swab platform suitable for automation, gram stains, traditional culture, and molecular assays. The study has highlighted a number of potential issues with the Copan eSwab™ in a direct comparison to the Σ-Transwab® Foam tip and the Σ-Transwab® PurFlock® tip.

Day	eSwab			Σ-Transwab® Foam			Σ-Transwab® PurFlock			Σ-Virocult®		
	Target Flu A	Average CT	Median CT	Target Flu A	Average CT	Median CT	Target Flu A	Average CT	Median CT	Target Flu A	Average CT	Median CT
1	inhibited	-	27.83	Detected	24.11	24.39	Detected	23.72	23.78	Detected	22.73	22.68
2	Detected	28.17		Detected	24.11		Detected	23.74		Detected	22.68	
3	Detected	27.64		Detected	24.69		Detected	23.89		Detected	22.59	
4	inhibited	-		Detected	24.57		Detected	23.79		Detected	22.60	
5	Detected	27.62		Detected	24.43		Detected	23.79		Detected	22.70	
Target Flu B	Average CT	Median CT	Target Flu B	Average CT	Median CT	Target Flu B	Average CT	Median CT	Target Flu B	Average CT	Median CT	
1	inhibited	-	29.70	Detected	25.12	24.82	Detected	24.20	24.07	Detected	23.01	23.01
2	Detected	29.67		Detected	24.65		Detected	23.88		Detected	23.00	
3	Detected	29.97		Detected	24.98		Detected	24.06		Detected	23.09	
4	inhibited	-		Detected	24.89		Detected	23.99		Detected	22.97	
5	Detected	29.55		Detected	24.66		Detected	24.31		Detected	22.93	
Target RSV	Average CT	Median CT	Target RSV	Average CT	Median CT	Target RSV	Average CT	Median CT	Target RSV	Average CT	Median CT	
1	inhibited	-	28.58	Detected	23.48	23.44	Detected	23.24	23.40	Detected	22.69	22.58
2	inhibited	-		Detected	22.94		Detected	23.38		Detected	22.46	
3	Detected	28.58		Detected	22.84		Detected	23.30		Detected	22.64	
4	Detected	-		Detected	23.73		Detected	23.48		Detected	22.44	
5	Detected	28.06		Detected	23.57		Detected	23.75		Detected	22.60	

Target	eSwab			Σ-Transwab® Foam			Σ-Transwab® PurFlock			Σ-Virocult®		
	Target	Average CT	Median CT	Target	Average CT	Median CT	Target	Average CT	Median CT	Target	Average CT	Median CT
Staphylococcus aureus	Detected	23.07	23.09	Detected	17.87	17.86	Detected	16.69	16.47	Detected	16.17	16.25
1	Detected	23.20		Detected	17.89		Detected	16.42		Detected	16.10	
2	Detected	23.13		Detected	17.93		Detected	16.43		Detected	16.26	
3	Detected	22.98		Detected	17.68		Detected	16.51		Detected	16.29	
4	Detected	23.08		Detected	18.02		Detected	16.56		Detected	16.30	
5	inhibited	-		Detected	17.56		Detected	16.91		Detected	15.44	15.32
Target Neisseria gonorrhoeae	Average CT	Median CT	Target Neisseria gonorrhoeae	Average CT	Median CT	Target Neisseria gonorrhoeae	Average CT	Median CT	Target Neisseria gonorrhoeae	Average CT	Median CT	
1	inhibited	-	18.52	Detected	17.15	17.44	Detected	16.96	16.94	Detected	15.53	
2	Detected	18.50		Detected	17.21		Detected	16.99		Detected	15.39	
3	Detected	18.71		Detected	17.49		Detected	17.01		Detected	15.34	
4	Detected	18.68		Detected	17.53		Detected	16.91		Detected	15.24	
5	inhibited	-		Detected	17.53		Detected	16.91		Detected	15.24	
Target Streptococcus pyogenes	Average CT	Median CT	Target Streptococcus pyogenes	Average CT	Median CT	Target Streptococcus pyogenes	Average CT	Median CT	Target Streptococcus pyogenes	Average CT	Median CT	
1	Detected	20.48	20.44	Detected	20.10	19.98	Detected	18.59	19.56	Detected	17.96	17.76
2	Detected	20.28		Detected	19.92		Detected	18.69		Detected	17.43	
3	Detected	20.43		Detected	19.90		Detected	18.89		Detected	17.90	
4	Detected	20.37		Detected	20.10		Detected	18.77		Detected	17.73	
5	Detected	20.90		Detected	19.49		Detected	18.56		Detected	17.64	

Khan, K. - An investigation of the compatibility of a new molecular transport medium (Σ-MM™) with a PCR analytical platform - Poster at the European Meeting on Molecular Diagnostics 11th Edition 2019

Table 1. CRE findings on Cepheid GeneXpert® PCR using Sigma Transwab®

Dilution	Samples tested for each organism	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours
10-1	2	IMP-1	Detected	Moderate Growth	Detected	Moderate Growth	Detected	Moderate Growth	Detected	Moderate Growth	Detected
10-2	2	Detected	Moderate Growth	Detected	Moderate Growth	Detected	Moderate Growth	Detected	Moderate Growth	Detected	Moderate Growth
10-3	2	Detected	Scanty Growth	Detected	Scanty Growth	Detected	Scanty Growth	Detected	Scanty Growth	Detected	Scanty Growth

Table 2. CRE findings on Cepheid GeneXpert® PCR using Sigma MM™

Dilution	Samples tested for each organism	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours	Target	Incubation culture growth results after 18 hours
10-1	2	IMP-1	Detected	No Growth	Detected	No Growth	Detected	No Growth	Detected	No Growth	No Growth
10-2	2	Detected	No Growth	Detected	No Growth	Detected	No Growth	Detected	No Growth	Detected	No Growth
10-3	2	Detected	No Growth	Detected	No Growth	Detected	No Growth	Detected	No Growth	Detected	No Growth

ZIEL

Demonstration of the compatibility of samples collected and inactivated with Σ-MM™ for testing on the Cepheid GeneXpert PCR Platform. The study was conducted in two phases: Phase 1: Σ-MM™ and Σ-Transwab® devices were spiked with specific dilutions of MRSA. Phase 2: Σ-MM™ and Σ-Transwab® devices were spiked with specific dilutions of CRE organisms (specific target sequences).

ERGEBNISSE

From the samples inoculated into Σ-MM™, nothing was recovered by culture, but all gave positive and correct results with GeneXpert (even low concentrations). From the samples inoculated into Σ-Transwab®, all were correctly recovered by culture, with the amount of growth appropriate to the concentration, and all gave positive and correct results with GeneXpert.

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Sample collected and stored in Σ-MM™ are compatible with GeneXpert. Σ-MM™ is effective at completely killing the sample, thus rendering it non-infectious. Σ-MM™ retains high quality bacterial DNA. Σ-Transwab® is also fully compatible with GeneXpert, and gave positive results even at low concentrations of target organisms.

Serobac Labordiagnostika Handelsgesellschaft m.b.H.
Neulinggasse 29 / Stiege 1 / Top 15, 1030 Wien
Tel.: +43 1 713 72 92, Fax: +43 1 712 37 99
E-Mail: serobac@serobac.at

www.serobac.at