



# Haemostasis Kits by Pentapharm



### **Pefakit® APC-R Factor V Leiden**

Pefakit® APC-R Factor V Leiden is a plasma based functional assay for the determination of resistance of Factor Va to inactivation by activated protein C (APC) caused by the factor V Leiden mutation (FV:Q<sup>506</sup>)

<b>Pack:</b>	3 x 2 ml	APC /RVV-V (+APC) Reagent
	3 x 2 ml	RVV-V (-APC) Reagent
	3 x 4 ml	PTA Reagent
	3 x 2 ml	Dilution Plasma
<b>REF:</b>	502-01	

### **Pefakit® APC-R Factor V Leiden Controls**

<b>Pack:</b>	3 x 1 ml	FV-L Negative Control
	3 x 1 ml	FV-L Heterozygous Control
<b>REF:</b>	502-21	

### **Pefakit® APC-R Factor V Leiden Homozygous Control**

<b>Pack:</b>	3 x 0.5 ml	FV-L Homozygous Control
<b>REF:</b>	502-23	

### **Pefakit® Reptilase® Time**

For the investigation of the last phase of blood coagulation. Due to its heparin insensitivity the Reptilase® Time can detect fibrinogen polymerization disorders even in the presence of heparin.

<b>Pack:</b>	3 x 1 ml	Reptilase® Time Reagent
<b>REF:</b>	800191	

### **Pefakit® in-TDT**

Determination of initial thrombin formation by intrinsic activation.

<b>Pack:</b>	2 x 2 ml	Reagent
	2 x 2 ml	Activator
	2 x 2,5 ml	Diluent
	1 x 1 ml	Calibrator
	1 x 1 ml	Control N
	1 x 1 ml	Control P
<b>REF:</b>	800236	

### **Pefakit® TAFI**

Pefakit® TAFI is a plasma based chromogenic assay for determination of **T**hrombin **A**ctivatable **F**ibrinolysis **I**nhibitor (TAFI) enzyme activity.

<b>Pack:</b>	2 x 4 ml	Activator
	2 x 4 ml	Start Reagent
	2 x 4 ml	Diluent
<b>REF:</b>	800186	

### **Pefakit® TAFI Calibrator and Controls**

<b>Pack:</b>	1 x 1 ml	Calibrator
	1 x 1 ml	Control 1
	1 x 1 ml	Control 2
<b>REF:</b>	800187	

### **Pefakit® PiCT® - Prothrombinase induced Clotting Time**

Pefakit® PiCT® is a plasma based functional assay for the determination of factor Xa and factor IIa inhibitors.

**Pack:** 3 x 2 ml PiCT® Activator  
3 x 2 ml PiCT® Start Reagent  
**REF:** 505-01

- **Applications for various haemostasis systems on request**
- **Applications for direct factor Xa inhibitors on request**

### **Pefakit® PiCT® Calibrators UFH**

**Pack:** 3 x 1 ml UFH Calibrator 1  
3 x 1 ml UFH Calibrator 2  
**REF:** 505-12

### **Pefakit® PiCT® Controls UFH**

**Pack:** 3 x 1 ml UFH Control 1  
3 x 1 ml UFH Control 2  
**REF:** 505-22

### **Pefakit® PiCT® Calibrators LMWH**

**Pack:** 3 x 1 ml LMWH Calibrator 1  
3 x 1 ml LMWH Calibrator 2  
**REF:** 505-11

### **Pefakit® PiCT® Controls LMWH**

**Pack:** 3 x 1 ml LMWH Control 1  
3 x 1 ml LMWH Control 2  
**REF:** 505-21

### **Pefakit® PiCT® Calibrators Hirudin**

**Pack:** 3 x 1 ml Hirudin Calibrator 1  
3 x 1 ml Hirudin Calibrator 2  
**REF:** 505-13

### **Pefakit® PiCT® Controls Hirudin**

**Pack:** 3 x 1 ml Hirudin Control 1  
3 x 1 ml Hirudin Control 2  
**REF:** 505-23